



Wiring diagram of the Control Panel - Connection of the taking-off sensors in control panel used for
MODEL PRO 500 (DC 310), MODEL PRO 700 (DC 350), DC 240, DC 290, DL 210, DL 270, DL 290, DL 320, DL 330, DL 340, DL 380,
MODEL PRO 600 (DA 270), DA 280, MODEL PRO 800 (DA 370), DA 380
 (1.4.2021)



BK - Black (čierna)
 BN - Brown (hnedá)
 BU - Blue (modrá)
 GN - Green (zelená)
 GNYE - Green/Yellow
 GY - Grey (sivá)
 OG - Orange (oranžová)
 PK - Pink (ružová)
 RD - Red (červená)
 VT - Violet (fialová)
 WH - White (biela)
 YE - Yellow (žltá)

Control panel with 5 instruments,
 lower hose delivery



Note: if dental unit not contain hygiene
 do not connect syringe switch

Control panel with 6 instruments,
 lower hose delivery



Note: if dental unit not contain hygiene
 do not connect syringe switch